


Neotech RAM Cannula[®]

Nasal Oxygen Cannula



FEATURES AND BENEFITS

- For use with low/high flow oxygen
- Can be used with or without humidity
- Soft, curved prongs
- Color coded sizes
- 15 mm oxygen tubing adapter included
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

<i>Cat No.</i>	<i>Item</i>	<i>Qty/Unit</i>
N4900	Micro Premie	10/box
N4901	Premie	10/box
N4902	Newborn	10/box
N4903	Infant	10/box
N4904	Small	10/box
N4905	Medium	10/box
N4906	Large	10/box
N4820	Oxygen Tubing Adapter	50/box



Comfortable for patient



Color coded sizes



Oxygen Adapter



#8,353,296; #D681,193; #9,027,562;
 Additional Patents Pending

Neotech RAM Cannula®

Nasal Oxygen Cannula

INDICATIONS FOR USE

The Neotech RAM Cannula is intended to be used in the administration of oxygen through the nares to pediatric (neonates, infants and children) patients.

DIRECTIONS FOR USE

Step 1

Select proper size cannula to fit patient.
Ideal prongs size will fill 60-80% of nares.

Note: Ensure that prongs do not fill nares completely.

Step 2

Attach 15 mm oxygen tubing adapter (provided) to RAM cannula. Attach cannula to oxygen source.

Note: Ensure all cannula connections are properly secured.

Step 3

Begin gas flow.

Step 4

Insert prongs into patient's nares.

Note: Allow small gap between patient's septum and base of prongs.

Step 5

Secure cannula to face with EZ-Hold® Tubing and Cannula Holder or comparable product.

Step 6

Replace Neotech RAM Cannula every 14 days or per hospital protocol, whichever is sooner.

GENERAL INFORMATION

Product Materials

Cannula tube: PVC

Cannula prongs: PVC

15mm adapter: PVC

O-ring: Silicone

Oxygen adapter: Acrylic

Made with natural rubber latex: No

Made with plasticizer DEHP: No

Size

N4900 Micro Preemie

Nasal prongs septal space: 2.5 mm

Nasal prongs outer diameter: 3 mm

N4901 Preemie

Nasal prongs septal space: 4.2 mm

Nasal prongs outer diameter: 3 mm

N4902 Newborn

Nasal prongs septal space: 4.8 mm

Nasal prongs outer diameter: 3.5 mm

N4903 Infant

Nasal prongs septal space: 5 mm

Nasal prongs outer diameter: 4 mm

N4904 Small

Nasal prongs septal space: 5.75 mm

Nasal prongs outer diameter: 5.25 mm

N4905 Medium

Nasal prongs septal space: 6.75 mm

Nasal prongs outer diameter: 5.75 mm

N4906 Large

Nasal prongs septal space: 7.75 mm

Nasal prongs outer diameter: 6.5 mm

Packaging

Individually packaged: Yes

Package type: Polyethylene bag

Dispenser box: Yes

Visit neotechproducts.com for the most up-to-date information.

Polyethylene bags are biodegradable. Boxes are printed with environmentally friendly soy ink on recyclable material.

CONTRAINDICATIONS, CAUTIONS & WARNINGS

The Neotech RAM Cannula is contraindicated for patients with nasal atresia or patients with facial structure deformities that prohibit adequate respiratory support.

Federal law restricts this device to sale by or on the order of a physician.

Improper selection of size, improper positioning or improper use may result in septal trauma or necrosis.

Always begin gas flow prior to inserting prongs into patient's nares.

Frequent observation of prongs position in patient's nares may be necessary

Cannula tubing can pose a potential strangulation hazard.

Do not use if package is damaged or open.

Do not soak, wash or sterilize.

Discontinue immediately if skin irritation occurs.

Single patient use only.



For more information or to order your FREE sample call
1-800-966-0500 or visit neotechproducts.com

Neotech Products LLC • 28430 Witherspoon Pkwy • Valencia, CA 91355, USA