Neotech RAM Cannula® Nasal Oxygen Cannula

REF Catalog No.

N4900 Micro Preemie
N4901 Preemie
N4902 Newborn
N4903 Infant
N4904 Small
N4905 Medium
N4906 Large

DIRECTIONS FOR USE
INSTRUCTIONS D'UTILISATION
NÁVOD K POUŽITÍ
GEBRAUCHSANWEISUNG
ISTRUZIONI PER L'USO
INSTRUCCIONES DE USO
ZALECENIA DOTYCZĄCE
UŻYTKOWANIA

 Manufactured by
NEOTECH PRODUCTS LLC
28430 Witherspoon Parkway, Valencia, CA 91355, USA
1-800-966-0500 • neotechproducts.com

EC **REP** EMERGO EUROPE
Prinsessegracht 20
2514 AP, The Hague, The Netherlands

Made in USA

D462 Rev E

English
Neotech RAM Cannula®
Nasal Oxygen Cannula

FEATURES & BENEFITS

- For use with low/high flow oxygen
- Can be used with or without humidity
- Soft, curved prongs
- Color coded sizes
- 15 mm oxygen tubing adapter included
- Not made with natural rubber latex or plasticizer DEHP
- Individually packaged

INDICATIONS FOR USE

The Neotech RAM Cannula is intended to be used in the administration of oxygen through the nares to pediatric (neonates, infants and children) patients.

DIRECTIONS FOR USE

Step 1 Select proper size cannula to fit patient. Ideal prongs size will fill 60-80% of nares.

Note: Ensure that prongs do not fill nares completely.

Step 2 Attach 15 mm oxygen tubing adapter (provided) to RAM cannula. Attach cannula to oxygen source.

Note: Ensure all cannula connections are properly secured.

Step 3 Begin gas flow.

Step 4 Insert prongs into patient's nares.

Note: Allow small gap between patient's septum and base of prongs.

Step 5 Secure cannula to face with EZ-Hold® Tubing and Cannula Holder or comparable product.

Step 6 Replace Neotech RAM Cannula every 14 days or per hospital protocol, whichever is sooner.

CONTRAINDICATIONS, CAUTIONS & WARNINGS

The Neotech RAM Cannula is contraindicated for patients with nasal atresia or patients with facial structure deformities that prohibit adequate respiratory support.

Federal law restricts this device to sale by or on the order of a physician.

Improper selection of size, improper positioning or improper use may result in septal trauma or necrosis.

Always begin gas flow prior to inserting prongs into patient's nares.

Frequent observation of prongs position in patient's nares may be necessary

Cannula tubing can pose a potential strangulation hazard.

Do not use if package is damaged or open.

Do not soak, wash or sterilize.

Discontinue immediately if skin irritation occurs.

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