NEOTECH®

Products LLC

Neotech RAM Cannula® Nasal Oxygen Cannula


REF

Catalog No.

N4900	N4904
N4901	N4905
N4902	N4906
N4903	

Manufactured by NEOTECH PRODUCTS LLC 28430 Witherspoon Parkway , Valencia, CA 91355, USA 1-800-966-0500 • neotechproducts.com

EMERGO EUROPE
Westervoortsejdijk 60
6827 AT, Arnhem, The Netherlands

Flow Resistance
N4900, N4901, N4902, N4903
0.446 cmH2O / L / min @ 2.5 LPM
N4904, N4905, N4906
0.230 cmH2O / L / min @ 15 LPM

English
Neotech RAM Cannula® Nasal Oxygen Cannula INDICATIONS FOR USE
The Neotech RAM Cannula is intended to be used for preterm and term neonates, infants, and pediatric patients who require supplemental oxygen, respiratory support or assistance in breathing in an ambulatory, hospital or institutional environment. It is used to provide oxygen with the following methods:
1. Heated or unheated, humidified low flow oxygen.
2. Heated, humidified high flow oxygen.
DIRECTIONS FOR USE
Step 1 Select proper size cannula by using RAM size guide provided in package. Ideal prongs size will fill approximately 80% of nares. Ensure that prongs do not fill nares completely.
Step 2 Attach 15 mm oxygen tubing adapter (provided) to RAM cannula. Attach cannula to oxygen source. Ensure all cannula connections are properly secured.
Step 3 Open and test the level of oxygen flow before putting the device on the patient.
Step 4 Insert prongs into patient’s nares and allow a small gap between patient’s septum and base of prongs.
Step 5 Secure cannula to face with EZ-Hold® Tubing and Cannula Holder or comparable product.
Step 6 Replace Neotech RAM Cannula every 14 days or per hospital protocol, whichever is sooner. Dispose of used RAM Cannula per hospital protocol.
CONTRAINdicATIONS, CAUTIONS & WARNINGS
The Neotech RAM Cannula is contraindicated for patients with nasal atresia or patients with facial structure deformities that prohibit adequate respiratory support.
Federal law restricts this device to sale by or on the order of a physician.
This device should only be used while the patient is under the continuous, direct supervision of healthcare professional and frequent observation of prongs position in patient’s nares is necessary

Improper selection of size, improper positioning or improper use may result in septal trauma or necrosis.
Cannula tubing can pose a potential strangulation hazard. Do not use the O-ring to secure the device to the patient.
Always begin gas flow prior to inserting prongs into patient’s nares. Do not use if package is damaged or open.
Do not soak, wash or sterilize.
Reuse of this device could result in cross contamination.
Discontinue immediately if skin irritation occurs.
Single patient use only.
European Union: MDR EU 2017/745
Any serious incident that has occurred in relation to this device should be reported to Neotech Products and the competent authority of the Member State in which the user and/or patient is established.

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