NEOTECH® Products LLC	
Part #	D492
Revision	M
Description	NeoSeal DFU
Artwork #	AW4492
Material # (PM, DM, LM) & description	DM466
Effective date	
CO #	

NEOTECH® Products LLC	
NeoSeal® Nasal Prong Seal	
REF	Catalog No.
	N420 Micro
	N421 Small
	N422 Large
	Manufactured by NEOTECH PRODUCTS LLC 28430 Witherspoon Parkway , Valencia, CA 91355, USA 1-800-966-0500 • neotechproducts.com
	EMERGO EUROPE Westervoortsejdkij 60 6827 AT, Arnhem, The Netherlands
	Medical Device
	
	
	D492 Rev M 2025.02.10

English
NeoSeal® Nasal Prongs Seal
INDICATIONS FOR USE
The NeoSeal is intended to create a seal around the nares and protect the skin of the nasal septum and nares. It is intended for use on pediatric (neonates, infants and children) patients.
DIRECTIONS FOR USE
Step 1 Peel NeoSeal from plastic liner.
Step 2 Slip over nasal prongs with adhesive facing cannula.
Note: Soft, smooth side without adhesive should contact patient's skin.
Step 3 Insert prongs per hospital protocol. Document placement of the NeoSeal in the patient’s medical record.
Note: To minimize risk of aspiration/ingestion, ensure the nasal cannula prongs and/or NeoSeal do not enter the infant’s mouth.
Step 4 Replace NeoSeal daily or per hospital protocol, whichever is sooner.
Note: Replace immediately if soiled or wet.
Step 5 Document removal and/or replacement in the patient’s medical record.
CAUTIONS & WARNINGS
Federal law restricts this device to sale by or on the order of a physician.
Choking hazard. NeoSeal should only be used while the patient is under the continuous, direct supervision of trained healthcare professionals.
Risk of ingestion or aspiration. When in use, ensure the nasal cannula prongs and/or NeoSeal do not enter the infant’s mouth. Replace immediately if soiled or wet.
Discontinue use immediately if skin irritation occurs.
Single patient use only.
European Union: MDR EU 2017/745
Any serious incident that has occurred in relation to this device should be reported to Neotech Products and the competent authority of the Member State in which the user and/or patient is established.

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