

# Making Products that Make a Difference



**Our mission and our motto are the same...  
*...Make a difference!***

Our unique products **make a difference** in the lives of babies, parents and healthcare practitioners throughout the world!

Our employees **make a difference** by providing friendly, prompt service and doing whatever it takes to provide 100% Customer Satisfaction!

Our quality and avoidance of back orders **make a difference** to the hospital or distributor in need of product!

We will continue to **make a difference** by introducing several new products a year!



RAM Cannula . . . . .	1	EZ-Hold . . . . .	13
Little Sucker . . . . .	2	NeoBridge . . . . .	14
Little Sucker Cover . . . . .	3	Hold-A-Line . . . . .	15
Suction Caddy . . . . .	4	NeoFlex . . . . .	16
Meconium Aspirator . . . . .	5	NeoBond . . . . .	17
NeoLead . . . . .	6	Lead Saver . . . . .	17
Micro NeoLead . . . . .	7	NeoSmile . . . . .	17
NeoBar for Neonates . . . . .	8	NeoGrip . . . . .	18
NeoBar for Pediatrics . . . . .	9	NeoPulse . . . . .	19
NeoShades with Tabs . . . . .	10	EZCare . . . . .	20
NeoShades with Headstrap . . . . .	11	EZCare Softouch . . . . .	21
Binasal Airway . . . . .	12		

**ALL PRODUCTS ARE LATEX, BPA AND PHTHALATE (DEHP) FREE**



# Neotech RAM Cannula™ **NEW**

A Simple, Revolutionary Interface for Nasal Respiratory Support



## FEATURES & BENEFITS

- Can deliver CPAP, PPV, IMV, and continuous oxygen
- Unique nasal cannula with universal 15mm adapter
- Soft, curved prongs currently available in three sizes
- Shorter, kink-resistant tubing allows for better pressures, better flows and less dead space
- Compatible with HFOV, T-Piece resuscitators, traditional wall oxygen, most ventilators and Bubble CPAP
- New oxygen adapter allows transition from NCPAP or NIPPV to high- or low-flow oxygen without having to change the cannula

Cat No.	Prod/Size	Patient Range	Qty/Unit
<b>Non Sterile</b>			
N4901NS	Preemie	< 1000 gm	10/box
N4902NS	Newborn	1000 - 2500 gm	10/box
N4903NS	Infant	> 2500 grams	10/box
N4820	Oxygen Adapter		50/box

[www.ramcannula.com](http://www.ramcannula.com)

Clinically proven as a safe, effective and gentle way to deliver NCPAP or NIPPV to babies.




Patents Pending **Made in USA**

Contact NEOTECH today: **1-800-966-0500** | **1-800-966-0585** FAX

# Little Sucker®

Oral-Nasal Suction Devices for Gentle Yet Effective Suctioning



## FEATURES & BENEFITS

- Single-handed suctioning
- Thumb port for intermittent suctioning
- Soft, flexible tip like a bulb syringe
- Use in NICU, PICU, Respiratory, Labor and Delivery
- Improve patient care and save nursing time
- Replace up to 3 products with 1 Little Sucker



Little Sucker with thumb port occluded.

Cat No.	Item	Qty/Unit
N204	Premie Little Sucker	50/box
N205	Standard Little Sucker	50/box
N224	Premie Nasal Tip	50/box
N225	Standard Nasal Tip	50/box



Patent #6,958,050; #7,635,361;  
#D474,269; #D590,056;  
Other International &  
Domestic Patents Pending

# Little Sucker® Cover **NEW**

Protective Cover for the Little Sucker



Little Sucker Cover with convenient cap off for thumb port.

## FEATURES & BENEFITS

- A safe and easy way to keep the Little Sucker clean between use
- May reduce risk of contamination
- Color helps you to quickly locate Little Sucker
- Shut off suction to reduce noise and save energy
- Easy to remove for immediate use

<i>Cat No.</i>	<i>Item</i>	<i>Qty/Unit</i>
N204C	Preemie w/Cover	50/box
N205C	Standard w/Cover	50/box
N215	Cover Only (Sucker not included)	50/box



Patents Pending **Made in USA**

All products are **Latex, BPA and Phthalate (DEHP) Free**

# Suction Caddy™

## Suction Device Holder to Store and Reuse the Neotech Little Sucker and Other Devices

### FEATURES & BENEFITS

- Attaches to most smooth surfaces; large suction cup provides strong hold and tab allows for easy removal
- Low price and easy to use
- Keep your suction device clean and close at hand
- Prevent tubing and devices from being misplaced or falling on the floor
- Cap off suction to reduce noise
- Write-on surface for patient information
- Ideal for holding the Neotech Little Sucker and other devices

Cat No.	Item	Qty/Unit
N245	Suction Caddy	40/box

Little Sucker placed in Suction Caddy and suction capped off.  
Little Sucker not included.




Made in USA Patent #7,601,001

# Neotech Meconium Aspirator<sup>®</sup>

Meconium Removal Device that Allows Intermittent Suctioning




THE FIRST  
& THE BEST



## FEATURES & BENEFITS

- Eliminates mouth suctioning of meconium
- Attaches to any endotracheal tube
- Thumb port allows for intermittent suctioning
- Featured in the Textbook on Neonatal Resuscitation
- Used worldwide in over 40 countries
- Available in yellow tint or original clear

Meconium Aspirator connected to ET tube and suction line with thumb port occluded.



Cat No.	Color	Qty/Unit
N0101	Clear	40/box
N0102	Yellow	40/box



All products are proudly Made in the USA

# NeoLead®

## Unique Electrodes with Hydrogel Center and Hydrocolloid Border for Long Lasting Leads

### FEATURES & BENEFITS

- All leads are pre-wired

#### All leads have hydrogel center

- Produces excellent tracings
- Safe on neonatal skin
- Won't dry out prematurely

#### All leads have hydrocolloid border

- Provides long term attachment
- Prevents gel from oozing out the edges
- Safe on neonatal skin
- Improves sensor stability




#### Radiolucent NeoLead (N305)

- Doesn't show up on x-rays

#### PedLead (N315)

- Two channel continuous monitoring
- Radiolucent materials allow for x-ray

Cat No.	Item	Size	Qty/Unit
N301	Standard	1" Circle	20 sets/box
N305	Radiolucent	1" Circle	20 sets/box
N315	Pediatric	1" Circle	20 sets/box



Made in USA Patent #6,795,722; Others Pending

# Micro NeoLead<sup>®</sup>

Unique Electrode with a Hydrogel Center and Hydrocolloid Ends



## FEATURES & BENEFITS

- Pre-wired
- No cumbersome band
- Safely secures to limbs without a strap that could reduce circulation
- Gentle hydrocolloid ends and hydrogel adhesives
- Long lasting and easy to use, even under high humidity
- Radiolucent materials allow for x-ray wherever they are positioned
- Safe on neonatal skin

<i>Cat No.</i>	<i>Item</i>	<i>Size</i>	<i>Qty/Unit</i>
N306	Micro	1" x 3/8"	20 sets/box



Micro NeoLeads applied to baby.



Patent #6,795,722; Others Pending

Made in USA

Contact NEOTECH today: 1-800-966-0500 | 1-800-966-0585 FAX

# NeoBar® for Neonates

## ET Tube Holder with Gentle Hydrocolloid Adhesive Tabs

### FEATURES & BENEFITS

- Attaches to patients with NeoBond® hydrocolloid
- Reduces extubations and helps prevent palate trauma
- Available in a variety of color coded sizes
- No tape near nose or mouth allows for better oral care
- Standardizes application procedures and prevents frequent tape changes
- Easily removed in case of emergency



<i>Cat No.</i>	<i>Size</i>	<i>Color</i>	<i>Qty/Unit</i>
N709	Mini	Purple	5/box
N710	Ultra	Yellow	5/box
N711	Micro	White	5/box
N712	Small	Green	5/box
N713	Large	Blue	5/box
N714	XL	Rose	5/box

Two box min. purchase. Ask about our assorted packs.

NeoBar attached to patient with gentle hydrocolloid tabs. ET tube secured to stabilization platform.



Patent #D434,496;  
#6,050,263; Others Pending

# NeoBar<sup>®</sup> for Pediatrics

ET Tube Holder with Acrylic or Gentle Hydrocolloid Adhesive Tabs



## FEATURES & BENEFITS

- Choice of NeoBond<sup>®</sup> hydrocolloid or medical grade acrylic adhesive tabs
- Reduces extubations and helps prevent palate trauma
- Available in a variety of color coded sizes
- No tape near nose or mouth allows for better oral care
- Standardizes application procedures and prevents frequent tape changes
- Easily removed in case of emergency

Cat No.	Size	Color	Qty/Unit
<b>With Medical Grade Acrylic Adhesive</b>			
N715F	Macro	Peach	5/box
N716F	Jumbo	Carbon	5/box
N717F	Mega	Brown	5/box
<b>With Gentle Hydrocolloid Adhesive</b>			
N715H	Macro	Clear	5/box
N716H	Jumbo	Gray	5/box
N717H	Mega	Tan	5/box

Two box min. purchase. Ask about our assorted packs.



Patent #D434,496;  
#6,050,263; Others Pending



All products are Latex, BPA and Phthalate (DEHP) Free

# NeoShades® with Tabs

Phototherapy Eye Shields with Cute Sunglass Design



## FEATURES & BENEFITS

- New, thicker, softer eye shield material with ocular pockets for superior patient comfort and safety
- Superior UV light blocking
- Improved tabs with ultra skin-friendly NeoBond® hydrocolloid
- Stronger, softer low profile hook



Cat No.	Size	Qty/Unit
N724	Micro	50/box
N725	Small	50/box
N726	Large	50/box

NeoShades with cute sunglasses design secured to patient with hydrocolloid tabs.



Patent #D410,021; Other International & Domestic Patents Pending

# NeoShades® with Headstrap

Phototherapy Eye Shields with NeoFoam® Headstrap



## FEATURES & BENEFITS

- New, thicker, softer eye shield material with ocular pockets for superior patient comfort and safety
- Superior UV light blocking
- Provides enhanced safety and less need for readjustment, reducing nursing time
- Patented NeoFoam keeps NeoShades in place and prevents slipping on hair and head molding

Cat No.	Size	Qty/Unit
N727	Micro	50/box
N728	Small	50/box
N729	Large	50/box



Special shape provides ideal fit for patients of all sizes.



Patent #D410,021; Other International & Domestic Patents Pending



All products are proudly Made in the USA

# Binasal Airway™

Silicone Nasopharyngeal NCPAP Prongs to Reduce Trauma

## FEATURES & BENEFITS

- Easy to insert and remove
- Made from soft, latex-free silicone
- Reduce septal and nares trauma
- Available in five different lengths and four different diameters
- Ideal for NCPAP and IMV
- Stays securely in place



**NO PIG NOSES**

Cat No.	OD MM (Fr)	Length CM	Qty/Unit
<b>Preemie</b>			
N82520	2.5 (7.5)	2	10/box
N83020	3.0 (9.0)	2	10/box
N82525	2.5 (7.5)	2.5	10/box
N83025	3.0 (9.0)	2.5	10/box
N82530	2.5 (7.5)	3	10/box
N83030	3.0 (9.0)	3	10/box
<b>Regular</b>			
N8254	2.5 (7.5)	4	10/box
N8259	2.5 (7.5)	9	10/box
N8304	3.0 (9.0)	4	10/box
N8309	3.0 (9.0)	9	10/box
N8354	3.5 (10.5)	4	10/box
N8359	3.5 (10.5)	9	10/box
N8404	4.0 (12)	4	10/box
N8409	4.0 (12)	9	10/box




# EZ-Hold®

## Tubing and Cannula Holder with Gentle Hydrocolloid Base and Secure Positioning Flaps



N756



N757

### FEATURES & BENEFITS

- Hydrocolloid base is gentle on delicate skin
- Fold over flap reduces stress on tubes
- Ideal for use with nasal cannulas and feeding tube
- Clear flap allows visualization of markings
- Easy to use
- Open and close repeatedly

<i>Cat No.</i>	<i>Size</i>	<i>Qty/Unit</i>
N756	Mini	100/box
N757	Standard	100/box




Mini EZ-Hold attached to patient with gentle hydrocolloid base. Cannula is secured by positioning flap.



Patent #D506,547 **Made in USA**

# NeoBridge®

## Umbilical Catheter Holder with Secure Positioning Flaps



### FEATURES & BENEFITS

- Attaches to patients with gentle NeoBond® Hydrocolloid
- Helps prevent dislodgement and standardizes procedure
- Available in three sizes
- Flaps can be opened and closed repeatedly
- Accommodates two catheters



NeoBridge secured to patient with hydrocolloid base while holding both UAC & UVC inside positioning flaps.

<i>Cat No.</i>	<i>Size</i>	<i>Approx Wt</i>	<i>Qty/Unit</i>
N700	Micro	800 or less	20/box
N701	Small	800-1200 gm	20/box
N702	Large	1200-2200 gm	20/box




Made in USA

Patent #5,370,627

# Hold-A-Line

## Catheter and Tubing Holder with Gentle Hydrocolloid Base and Secure Positioning Flap



### FEATURES & BENEFITS

- Easy to open and close repeatedly
- Holds catheters, IVs and various tubes
- Hydrocolloid base is gentle on skin
- Unique shape allows for larger umbilical stumps
- Clear flap allows visualization of markings
- Standardizes securement procedures
- Relieves stress on tubes

<i>Cat No.</i>	<i>Size</i>	<i>Qty/Unit</i>
N705	Standard	40/box

Hold-A-Line securely holds line in place and helps to prevent tube dislodgement.



Patents Pending **Made in USA**

# NeoFlex™ **NEW**

## Silicone Gel Adhesive Roll

### FEATURES & BENEFITS

- Studies show NeoFlex is safe on fragile skin
- Flexible and easy to apply, adjust or remove
- Great as a skin barrier
- Use to secure limb boards, and various other uses
- Ideal for micro preemies and sensitive skin
- Conveniently sized rolls for single patient use



<i>Cat No.</i>	<i>Size</i>	<i>Qty/Unit</i>
N125	3/4" wide x 3' long	20 rolls/box



Arm secured to limb board with NeoFlex. Clear material allows you to monitor skin.



Made in USA Patents Pending

# NeoBond®

Hydrocolloid Strips Can be Used in Place of Harmful Tape



- Provides long term attachment but prevents epidermal stripping
- Leaves no residue on skin

<i>Cat No.</i>	<i>Size</i>	<i>Qty/Unit</i>
N741	Standard*	50 sets/box
*Each NeoBond has two 1/2" x 6" strips		

# Lead Saver™

Electrode or Temp Probe Cover Securing Device



- Prevents electrodes from sliding
- Skin-friendly hydrocolloid improves sensor stability
- Prevents gel from oozing out edges
- Also use on temp probe covers
- Works on circular or square leads

<i>Cat No.</i>	<i>Size</i>	<i>Qty/Unit</i>
N325	Standard	300/box

# NeoSmile™

Hydrocolloid Temperature Probe Cover



- Temperature probe cover that attaches to patients with NeoBond hydrocolloid
- Resists skin breakdown from moisture or heat
- Doesn't fall off prematurely

<i>Cat No.</i>	<i>Size</i>	<i>Qty/Unit</i>
N731	1" Dia.	150/box

# NeoGrip™

## Multi-Purpose Tubing and Cable Holders

### FEATURES & BENEFITS

- Leaves no residue on tubes, cables or wires
- Holds suction tubes, no bunching or twisting
- Attaches to patient for better mobility
- Use in any hospital department
- Keep your equipment with you at all times
- Great for use in NICU
- Ideal for use during transport
- N603 can be sterilized for use in OR

Cat No.	Color	Qty/Unit
N601	Teal (VersaGrip)	40/box
N602	Pink/Blue	40/box
N603	White (sterilizable)	40/box




Tubes maintained by NeoGrip and attached to mother. Ideal for Kangaroo Care.



# NeoPulse™

Lightweight, Adjustable and Disposable Pulse Oximeter Wrap



NeoPulse gently wrapped around foot.  
Pulse oximeter not included.

## FEATURES & BENEFITS

- Made with unique, patented NeoFoam®
- Soft and ultra skin-friendly
- Easily adjustable to fit any size
- Lightweight and comfortable
- Provides stable, secure attachment that won't slip
- Breathable material resists development of foot odor
- Wicks away perspiration

Cat No.	Size	Dimensions	Qty/Unit
N781	Small	3/4" x 3"	100/box
N782	Large	1" x 5"	100/box



Patent Pending **Made in USA**

# EZCare™

## Tracheostomy Tube Holder

### FEATURES & BENEFITS

- Neoprene material is soft, lightweight and comfortable, and won't irritate the skin
- Provides stable attachment and reduces chance of accidental dislodgement
- Moisture resistant material allows patient to be bathed without removing holder
- Stretchable material allows for cough reflex
- Colorful straps help resist staining



Cat No.	Strap Size	Neck Size	Qty/Unit
N905	5"	7-8"	5/pack
N906	6"	8-9"	5/pack
N907	7"	9-10"	5/pack
N908	8"	10-11"	5/pack
N909	9"	11-12"	5/pack
N910	10"	12-13"	5/pack
*Tracheostomy tube not included			



Pink  
(PK)



Camouflage  
(CO)



Blue  
(BL)



Hawaiian  
(HN)



Made in USA Patents Pending

# EZCare Softouch™

## Disposable Tracheostomy Tube Holder



### FEATURES & BENEFITS

- Soft, skin-friendly NeoFoam™ won't slip
- Disposable, perfect for hospitals
- Lightweight, comfortable and breathable
- Easily adjusts for any neck size

<i>Cat No.</i>	<i>Size</i>	<i>Qty/Unit</i>
N924	Small (9")	50/box
N925	Medium (13")	50/box
N926	Large (18")	50/box
<b>*Tracheostomy tube not included</b>		



Patents Pending **Made in USA**

All products are Latex, BPA and Phthalate (DEHP) Free

# NEOTECH™

## Products, Inc.

27822 Fremont Court • Valencia, CA 91355 USA



**1-800-966-0500**  
[www.neotechproducts.com](http://www.neotechproducts.com)

**All Products are Latex, BPA  
and Phthalate (DEHP) Free**



**Friend us on Facebook**  
Neotech Products



**Follow us on Twitter**  
Neotechproducts



**Please Recycle**



**eco**  
ecoalition.org

