

Product Information Little Sucker®

- ◆ **N204 Little Sucker® - Premie** **50/case**
- ◆ **N205 Little Sucker® - Standard** **50/case**
- ◆ **N225 Nasal Tip Little Sucker® - Nasal** **50/case**
- ◆ **N224 Premie Nasal Tip Little Sucker® - Nasal** **50/case**
 - ◆ Oral - Nasal Suction Devices for Gentle Yet Effective Suctioning

✓ **Features and Benefits:**

- ✓ Singled-handed Suctioning
- ✓ Thumb port of intermittent suctioning
- ✓ Soft, flexible tip like a bulb syringe
- ✓ Latex-free
- ✓ Improve patient care and save nursing time
- ✓ Replace up to 3 products with 1 Sucker
- ✓ Use in NICU, PICU, Respiratory, L&D, ER, EMS, Home Health

➤ **DFU**

- Connect barbed end to standard suction line
- Set suction pressure between 80-100mm Hg or per hospital protocol
- Occlude thumb port to begin suctioning
- Gently suction mouth or nose with intermittent suctioning
- Suction sterile water through device after each use. Store between use in Suction Caddy.
- **Replace every 24 hours**