

NEOTECH™

Products, Inc.

27822 Fremont Court • Valencia, CA 91355 USA



1-800-966-0500

www.neotechproducts.com

All Products are Latex, BPA
and Phthalate (DEHP) Free

 Friend us on Facebook
Neotech Products

 Follow us on Twitter
Neotechproducts



NEOTECH™



**Perhaps You
Should Cover Up**
Introducing the Little Sucker Cover

Call 1-800-966-0500 or visit www.neotechproducts.com

Little Sucker®

Oral-Nasal Suction Devices for Gentle Yet Effective Suctioning



N204



N205



N224



N225

**NEW
ONE-PIECE
NASAL TIP
DESIGN!**


FEATURES & BENEFITS


- Single-handed suctioning
- Thumb port for intermittent suctioning
- Soft, flexible tip like a bulb syringe
- Use in NICU, PICU, Respiratory, Labor and Delivery
- Improve patient care and save nursing time
- Replace up to 3 products with 1 Little Sucker

Cat No.	Item	Qty/Unit
N204	Premie Little Sucker	50/box
N205	Standard Little Sucker	50/box
N224	Premie Nasal Tip	50/box
N225	Standard Nasal Tip	50/box



Premie Nasal Tip connected to suction line with thumb port occluded.



Patent #6,958,050; #7,635,361; #D474,269;
#D590,056; Other International & Domestic
Patents Pending  **Made in USA**

Little Sucker® Cover **NEW**

Protective Cover for the Little Sucker

FEATURES & BENEFITS

- A safe and easy way to keep the Little Sucker clean between use
- May reduce risk of contamination
- Color helps you to quickly locate Little Sucker
- Shut off suction to reduce noise and save energy
- Easy to remove for immediate use

Cat No.	Item	Qty/Unit
N204C	Preemie With Cover	50/box
N205C	Standard With Cover	50/box
N215	Cover Only (Little Sucker not included)	50/box



Little Sucker Cover with convenient cap off for thumb port.



Made in USA Patents Pending

Suction Caddy™

Suction Device Holder to Store and Reuse the Neotech Little Sucker and Other Devices



Little Sucker placed in Suction Caddy and suction capped off. Little Sucker not included.

FEATURES & BENEFITS

- Attaches to most smooth surfaces; large suction cup provides strong hold and tab allows for easy removal
- Low price and easy to use
- Keep your suction device clean and close at hand
- Prevent tubing and devices from being misplaced or falling on the floor
- Cap off suction to reduce noise
- Write-on surface for patient information
- Ideal for holding the Neotech Little Sucker and other devices

Cat No.	Item	Qty/Unit
N245	Suction Caddy	40/box

