

## Strap-Tote™

Affordable Monitor Protection

### FEATURES & BENEFITS

- Quick and easy to use
- Economical choice for monitor protection
- Velcro closure helps protect monitor from foods and liquids
- 5" x 6.5" Pocket

<i>Cat No.</i>	<i>Dimensions</i>	<i>Qty/Unit</i>
A-506	5" x 6.5"	100/box



Made in USA

1-800-966-0500 | [www.neotechproducts.com](http://www.neotechproducts.com)

27822 Fremont Court, Valencia, CA 91355

# Strap-Tote™ Directions For Use

## **Step 1**

- Open poly bag, remove pouch, unfold strap.

## **Step 2**

- Insert monitor and connect lead wires to monitor and patient.

## **Step 3**

- Hang strap around neck of the patient.