


## Valu-Tote™

Affordable Monitor Protection

### FEATURES & BENEFITS

- Our most economical pouch
- Secure clip for easy attachment
- Cost effective monitor protection
- 4.75" x 8.75" Pouch

Cat No.	Dimensions	Qty/Unit
A-507	4.75" x 8.75"	100/box



1-800-966-0500 | [www.neotechproducts.com](http://www.neotechproducts.com)

27822 Fremont Court, Valencia, CA 91355

# Valu-Tote™ Directions For Use

## ***Step 1***

- Insert monitor into pouch.

## ***Step 2***

- Connect lead wires.

## ***Step 3***

- Seal zipper on pouch.

## ***Step 4***

- Secure clip to patient gown.