

Versa-Tote™

Affordable Monitor Protection

FEATURES & BENEFITS

- Water resistant
- Secure clip for easy attachment
- Flap protects monitor from foods and liquids
- 5" x 9" Pouch

<i>Cat No.</i>	<i>Dimensions</i>	<i>Qty/Unit</i>
A-504	5" x 9"	100/box



Made in USA

1-800-966-0500 | www.neotechproducts.com

27822 Fremont Court, Valencia, CA 91355

Versa-Tote™ Directions For Use

Step 1

- Open pouch and insert transmitter with lead wires connected.

Step 2

- Seal Velcro® closure on pouch

Step 3

- Clip on the patient's gown with flap.
(Velcro side facing patient.)