

## Telemetry Totes



### FEATURES AND BENEFITS

#### Neotech Clip-Tote

- Our most popular pouch
- Water resistant
- Secure clip for easy attachment
- Allows for patient mobility

#### Neotech Valu-Tote

- Our most economical pouch
- Secure clip for easy attachment
- Allows for patient mobility

Cat No.	Product	Dimensions	Qty/Unit
A501	Clip-Tote	5" x 8"	100/box
A507	Valu-Tote	4.75" x 8.75"	100/box



Made in USA

# Telemetry Totes

## DIRECTIONS FOR USE

---

### Step 1

Open tote at zipper.

### Step 2

Insert telemetry transmitter into tote.

### Step 3

Thread lead wires through slit and connect to transmitter.

### Step 4

Seal zipper.

### Step 5

Secure clip to patient gown or robe.



For more information or to order your FREE sample call  
**1-800-966-0500** or visit [neotechproducts.com](http://neotechproducts.com)

Neotech Products LLC • 28430 Witherspoon Pkwy, Valencia, CA 91355